

Development and Evaluation of a Label Friendly Film Coating System Containing Naturally Derived Opacifier

Authors: Tara Moffa, Sydney Badger, Tom Mehaffey, Manish Ghimire, Ali Rajabi-Siahboomi

Colorcon, Inc. Harleysville, PA 19438, USA | AAPS Poster Reprint 2022

Introduction

Opacifiers are included in film coatings to provide opaqueness, help color matching, and protect active ingredients from light-induced degradation. Titanium dioxide (TiO₂), classified as E171 when used as a pigment in foods, has been the preferred opacifier in the industry. However, the use of TiO₂ as a food additive has been suspended in Europe, requiring companies to explore alternative materials. Mineral-based opacifiers such as calcium carbonate or magnesium carbonate have lower opacity than TiO₂ and are currently being used as an alternative to TiO₂. Formulators have achieved high opacity through optimization of the film coating within these mineral-containing systems¹. Our recent studies have shown that the alkaline nature of minerals can affect the color of naturally derived pigments. In addition, they may interact with some active ingredients. The current study aimed to evaluate the opacity of a novel fully formulated Nutrafinish[®] Label Friendly System using natural, plant-based opacifiers.

Methods

Nutrafinish[®] immediate release film coatings are developed specifically for nutritional and dietary supplement products and the formulations shown in Table 1 were used in this study. Nutrafinish high opacity and label friendly formulations were optimized to deliver high levels of opacity using calcium carbonate (CaCO₃) and with a natural opacifier, respectively. Their performance will be compared to a formulation containing TiO₂.

Table 1. Evaluated Nutrafinish Formulations

Nutrafinish System	Main Polymer	Opacifiers
Nutrafinish TiO ₂	HPMC	TiO ₂
Nutrafinish High Opacity	HPMC	CaCO ₃
Nutrafinish Label Friendly	HPMC	Plant based
Experimental	HPMC	Plant based

The formulations were coated onto 3.0 kg of multivitamin tablets (1000 mg) using a fully perforated 15" side-vented coating pan (Labcoat II, O'Hara Technologies Inc., Canada) to a 5% weight gain (WG) with samples taken at 3%, 4% and 5% WG to evaluate opacity. In addition to these systems, tablets were coated using anthocyanins in combination with either a calcium carbonate or a plant-based excipient as an opacifier.

Table 2. Coating Parameters

Parameter	Nutrafinish TiO ₂ / Experimental	Nutrafinish High Opacity/Label Friendly
Batch Size (g)	3000	3000
Spray Rate (g/min)	18 - 20	18-20
Dispersion Solids Content (%w/w)	20	20
Bed Temperature (°C)	45	45
Inlet Air Temperature (°C)	58	58
Air Flow (cfm/m ³ /hr)	165	165
Number of Spray Guns	1	1
Pan speed (rpm)	8 - 10	8 - 10
Atomization air (psi)	20	20
Pattern air (psi)	20	20

Results

Figure 1 shows uncoated and coated multivitamin tablets at 3, 4, and 5% coating weight gain.

Nutrafinish with TiO₂ showed the best opacity with uniform coverage at 3% WG (Figure 1b).

The best TiO₂ free option (Nutrafinish High Opacity) uses CaCO₃ as the opacifier, and uniformity can be achieved at 3% WG (Figure 1c).

When TiO₂ was directly replaced with a plant-based opacifier, the opacity was significantly decreased, and the coating was unable to hide the spotting of multivitamin tablets even at 5% WG (Figure 1d).

The novel optimized Nutrafinish Label Friendly formulation using a natural opacifier shows good opacity that can be achieved at 3-4% WG, without the use of calcium carbonate (Figure 1e).

Figure 1. (a) Uncoated Multivitamin, (b) Multivitamin Tablets coated with Nutrafinish TiO₂, (c) Nutrafinish High Opacity, (d) Experimental, and (e) Nutrafinish Label Friendly



Figure 2. Multivitamins Coated with Nutrafinish Containing Anthocyanin Pigments



Nutrafinish containing calcium carbonate has a high pH of ~ 8.5 causing the film coating color to turn brown. Nutrafinish Label Friendly has a neutral, pH of ~ 7 which allows it to develop a vibrant red color.

This demonstrates the advantage associated with calcium carbonate-free opacifiers for the color development of film coating systems with natural pigments.

Figure 3. Multivitamins Coated with Nutrafinish Label Friendly Systems Containing a Plant Based Opacifier



Conclusions

A novel Nutrafinish Label Friendly coating was developed using a naturally derived plant-based opacifier that replaced commonly used mineral components such as TiO₂ or calcium carbonate.

This system was designed specifically for the nutritional and dietary supplement market to meet consumer demand for a clean-label choice.

The advantage of a naturally derived opacifier within coating systems using non-synthetic pigments such as anthocyanins was shown to result in enhanced and vivid coloring compared to mineral-based formulations.

References

1. Badger S, Ghimire, M and Rajabi-Siahboomi, A. Impact of Opacifier Type in a Film Coating Formulation on Photostability of Tablet Ingredients. AAPS 360 (2021). T1430-07-37.

The information contained herein, to the best of Colorcon, Inc.'s knowledge is true and accurate. Any recommendations or suggestions of Colorcon, Inc. with regard to the products provided by Colorcon, Inc. are made without warranty, either implied or expressed, because of the variations in methods, conditions and equipment which may be used in commercially processing the products, and no such warranties are made for the suitability of the products for any applications that you may have disclosed. Colorcon, Inc. shall not be liable for loss of profit or for incidental, special or consequential loss or damages.

Colorcon, Inc. makes no warranty, either expressed or implied, that the use of the products provided by Colorcon, Inc., will not infringe any trademark, trade name, copyright, patent or other rights held by any third person or entity when used in the customer's application.

For more information, contact your Colorcon representative or call:

North America
+1-215-699-7733

Europe/Middle East/Africa
+44-(0)-1322-293000

Latin America
+54-11-5556-7700

India
+91-832-6727373

China
+86-21-61982300

You can also visit our website at www.colorcon.com