

Novel Authentication Technology Using Molecular Tags as a PCID in Solid Oral Dosage Forms

Authors: B. Prusak¹, D. To¹, L. Jung¹, M. Hogan² and A. Rajabi-Siahboomi¹

¹Colorcon, Inc. Harleysville, PA, 19438 USA

²Applied DNA Sciences, Inc., Stony Brook, NY, USA | AAPS Poster Reprint 2019

Introduction

A review, by The World Health Organization (WHO), of publications between 2007-2016 estimates that 10.5% of medicines are substandard or falsified and could account for loss of \$30.5 billion of pharmaceutical sales in low- and middle-income countries.¹ Medicines impacted cover a wide range of treatment categories including oncology, contraceptives, antibiotics, vaccines and other life-saving medicines.² The FDA has issued guidance to address this issue by incorporating physical or chemical identifiers (PCID) into solid dosage forms.³ Positive detection of the PCID can provide authentication and traceability to individual dosage forms. In this study, an Opadry[®] complete film coating system was used to carry a DNA-based molecular tag (SigNature[®], Applied DNA Sciences, Inc., USA), as a PCID covert authentication platform. This DNA tag is considered as a “molecular bar code”, enabling identification to the source, which could be a product type, or other meaningful attributes.

Methods

SigNature molecular tag was incorporated into samples of blue, gray and white pigmented Opadry film coating, then compared with untagged samples. Opadry powder samples were tested for the presence of the tag with SigNify[®] Reagent Mix and a SigNify[®] IF portable reader (Applied DNA Sciences, Inc., NY, USA) using real-time polymerase chain reaction (rt-PCR).

Tagged and untagged film coating samples were coated onto placebo tablets (10mm, bi-convex round), to a 3% weight gain using a Labcoat I (O’Hara Technologies, Inc., Ontario, Canada) with the processing conditions described in **Table 1**.

Samples of the film coated tablets were prepared and tested for the presence of the DNA tag, using the rt-PCR, as described above. Tablet appearance was compared analytically by measuring color difference using a DataColor600 reflection spectrophotometer (DataColor, Inc., USA), surface gloss using a Model 803A Surface

Analysis System (TRICOR Systems, Inc., USA) and surface roughness using a PS50 Optical Profilometer (Nanovea, Inc., USA). Limit of CIELAB total color difference (DE) were defined as 2.5, 2.0 and 1.5 for blue, gray and white samples, respectively. The performance of both tagged Opadry dry powder and film coated tablets under accelerated storage condition (40°C/75% RH) were also evaluated.

Table 1. Opadry Coating Process

Parameter	Opadry
Charge (kg)	1
Solids content (%w/w)	20
Spray rate (g/min)	8
Bed temp (°C)	45
Inlet air temperature	64
Air flow (cfm/m ³ /hr)	125/212
No. of guns	1
Gun type	VAU
Pan speed (rpm)	20
Atomization air (psi/bar)	20/1.4
Pattern air (psi/bar)	20/1.4

Results


Opadry samples were analyzed by rt-PCR and the resulting cycle threshold (Cq) was compared against negative and positive controls, as shown in **Figure 1a**. A higher Cq value indicates more cycle time needed to amplify DNA tag to threshold and therefore lower quantity of DNA present. In this study, any Cq value >25 were considered as non-detectable. The negative control and untagged Opadry samples, resulted in a Cq >30, indicating no detection of the tag. In comparison, the tagged Opadry samples and positive control had Cq values of 10-25 confirming presence of the tag.

Tablets coated with tagged and untagged Opadry were also tested by rt-PCR and indicated a similar trend. Untagged tablets and the negative control provided no detection (Cq values of greater than 30) while tagged tablets and positive controls resulted in reliable detection, as shown in **Figure 1b**. Comparisons of coated tablets using appearance

differentiation techniques: imaging, gloss, surface roughness, and color difference are shown in **Figures 2-4** and **Table 2**, respectively. These appearance techniques were not capable of identifying significant differences between tagged and untagged tablets.

Figure 1. Cycle Threshold (Cq) for (a) Opadry Powder and (b) Opadry Coated Tablets

(a) Opadry Powder



(b) Coated Tablets




Figure 2. Comparison of Images for Tablets Coated with Tagged and Untagged Opadry

Opadry Blue



Opadry Gray



Opadry White



Figure 3. Comparison of Gloss for Tablets Coated with Tagged and Untagged Opadry




Figure 4. Comparison of Surface Roughness of Tablets Coated with Tagged and Untagged Opadry





Table 2. Comparison of Color Difference of Tablets Coated with Tagged and Untagged Opadry


Parameter	Opadry Blue		Opadry Gray		Opadry White	
	Untagged	Tagged	Untagged	Tagged	Untagged	Tagged
Color Difference (DE)	Standard	0.16	Standard	0.43	Standard	0.47

Figure 5. Impact of Storage at 40°C/75% RH for 3 months on Cycle Threshold (Cq) for (a) Opadry Powder and (b) Opadry Coated Tablets

(a) Opadry Powder



(b) Coated Tablets



Conclusions

The SigNature molecular tag was successfully incorporated into an Opadry film coating as a covert PCID. This is an approach to uniquely identify tablets and capsules, and through authentication, protect against substandard and falsified medicines. The presence of tag could not be detected by appearance or common analytical methods; however, detection was confirmed by rt-PCR when used in conjunction with the matching SigNify Reagent Mix. Tagged samples of Opadry powder and coated tablets were differentiated from the untagged versions, indicating an excellent method for an on-tablet PCID authentication technology platform.

References

1. A Study on the Public Health and Socioeconomic Impact of Substandard and Falsified Medical Products. <https://www.who.int/medicines/regulation/ssffc/publications/Layout-SEstudy-WEB.pdf?ua=1> accessed April 2019
2. WHO Global Surveillance and Monitoring System for Substandard and Falsified Medical Products. https://www.who.int/medicines/regulation/ssffc/publications/GSMS_Report_layout.pdf?ua=1 accessed April 2019
3. Guidance for Industry - Incorporation of Physical-Chemical Identifiers into Solid Oral Dosage Form Drug Products for Anticounterfeiting. <https://www.fda.gov/downloads/drugs/guidances/ucm171575.pdf> accessed April 2019

The information contained herein, to the best of Colorcon, Inc.'s knowledge is true and accurate. Any recommendations or suggestions of Colorcon, Inc. with regard to the products provided by Colorcon, Inc. are made without warranty, either implied or expressed, because of the variations in methods, conditions and equipment which may be used in commercially processing the products, and no such warranties are made for the suitability of the products for any applications that you may have disclosed. Colorcon, Inc. shall not be liable for loss of profit or for incidental, special or consequential loss or damages.

Colorcon, Inc. makes no warranty, either expressed or implied, that the use of the products provided by Colorcon, Inc., will not infringe any trademark, trade name, copyright, patent or other rights held by any third person or entity when used in the customer's application.

You can also visit our website at www.colorcon.com

For more information, contact your Colorcon representative or call:

North America
+1-215-699-7733

Europe/Middle East/Africa
+44-(0)-1322-293000

Latin America
+54-11-5556-7700

India
+91-832-6727373

China
+86-21-61982300