



Wet Granulation of Acetaminophen with Starch 1500®

Technical Data

Objective

To demonstrate the combined binding and disintegration properties of Starch 1500 in a low shear wet granulation process. Acetaminophen powder was chosen as the example active due to its high dose, poor flow, and compaction properties. The use of Starch 1500 in the granulation was examined in two ways: As a dry mix with the acetaminophen using water as the granulation binder (Formula A), and with some of the Starch 1500 dispersed in the water to be used as the binder liquid (Formula B). The binding properties of Starch 1500 were also compared to PVP (Polyvidone), a commonly used wet granulation binder (Formula C).

Granulation Process

The granulations were conducted on a laboratory scale using an 8 qt. Hobart Planetary Mixer.

Formulation	А	В	С
Acetaminophen	85.10	85.10	85.10
Starch 1500 (dry)	14.65	11.73	9.65
Starch 1500 (in water)	-	2.92	-
PVP K 29/32 (in water)	-	-	5.00
Magnesium stearate	0.25	0.25	.25
Granulation Conditions	А	В	С
Hobart speed setting	1	1	1
Dry mix time (min.)	4	4	4
Wet mass time (min.)	5	5	5
Binder used	water alone	water alone + Starch 1500	water + PVP
Binder concentration (% solids)	-	20.0	18.2
Wet screening (mesh)	12	12	12
Drying conditions	Α	В	С
Glatt GPCG-3 Inlet air (deg.C.)	65	65	65
Drying time (min.)	21	27	35
Final product temperature (deg.C.)	40	40	40
Final granulation % L.O.D.	1.4	1.2	1.2
Blending - 8 qt. "V" blender	Α	В	С
Magnesium stearate blend time (min.)	3	3	3


Partially Pregelatinized Maize Starch

Compaction Process


Each of the granulations were compressed on a 10 station rotary tablet press using size B, 3/8" standard concave tooling, to a total tablet weight of 382.0 mg. Tablet samples were taken at 6 different compaction forces between 6 and 15 kN of force.

Compaction Profiles

The compaction profiles for the three granulations were very similar and each of the three formulations produced robust tablets. Adding some of the Starch 1500 to the water for granulation (Formula B) slightly increased the overall tablet hardness. The batch with PVP as the binder (Formula C) produced slightly harder tablets. However, upon disintegration testing in water the clear advantage of Starch 1500 as a binder was seen.



Disintegration Results



Conclusions

Starch 1500 exhibited dual functionality in this formulation. As a wet granulation binder it produced tablets with similar hardness to PVP. As a disintegrant it significantly outperformed PVP, which acutally caused a delay in disintegration.



World Headquarters

Colorcon

415 Moyer Blvd., P.O. Box 24, West Point, PA 19486-0024

Tel: 215-699-7733 Fax:215-661-2605 Website: www.colorcon.com e-mail: info@colorcon.com

Locations United States	Telephone	Facsimile	Locations Asia/Pacific	Telephone	Facsimile
Santa Ana, California	714-549-0631	714-549-4921	Singapore	65-6438-0318	65-6438-0178
Indianapolis, Indiana	317-545-6211	317-545-6218	Fuji-gun, Shizuoka, Japan	81-5-4465-2711	81-5-4465-2730
Humacao, Puerto Rico	787-852-3815	787-852-0030	Shanghai, China	86-21-5442-2222	86-21-5442-2229
			Goa, India	91-832-288-3434	91-832-288-3440
			Seoul, Korea	82-2-2057-2713	82-2-2057-2179
Europe			Latin America		
Dartford, Kent, England	44-1322-293000	44-1322-627200	Buenos Aires, Argentina	54-11-4552-1565	54-11-4552-3997
Dartford, Kent, England Bougival,, France	44-1322-293000 33-1-3082-1582	44-1322-627200 33-1-3082-7879	Buenos Aires, Argentina Cotia, Brasil	54-11-4552-1565 55-11-4612-4262	54-11-4552-3997 55-11-4612-3307
			, 3		
Bougival,, France	33-1-3082-1582	33-1-3082-7879	Cotia, Brasil	55-11-4612-4262	55-11-4612-3307
Bougival,, France Idstein, Germany	33-1-3082-1582 49-6126-9961-0	33-1-3082-7879 49-6126-9961-11	Cotia, Brasil Bogota, Colombia	55-11-4612-4262 571-418-1202	55-11-4612-3307 571-418-1257
Bougival,, France Idstein, Germany Gallarate, Italy	33-1-3082-1582 49-6126-9961-0 39-0331-776932	33-1-3082-7879 49-6126-9961-11 39-0331-776831	Cotia, Brasil Bogota, Colombia Caracas, Venezuela	55-11-4612-4262 571-418-1202 58-212-442-4819	55-11-4612-3307 571-418-1257 58-212-442-8724

The information contained herein, to the best of our knowledge is true and accurate. Any recommendations or suggestions are made without warranty or guarantee, since the conditions of use are beyond our control. Any information contained herein is intended as a recommendation for use of our products so as not to infringe on any patent.

